

A = Biopor[®] AB
40 Shore

B = Biopor[®] AB
70 Shore

C = Biopor[®] AB
fluoreszent
40 Shore

D = Biopor[®] AB
monitoring
40 Shore



00
farblos-transparent
A B



01
rötlich-transparent
A B



12
rot-opak
A B



14
orange-opak
A B



18
orange-fluoreszent
C



15
rosa-opak
A



16
rot-fluoreszent
C



17
pink-fluoreszent
C



81
himbeerrot-opak
B



82
erdbeerrot-opak
B



83
fuchsia-magenta-opak
B



19
lila-opak
A B



27
azurblau-opak
B



22
blau-opak
A B



26
blau-fluoreszent
C



32
gelb-opak
A B



34
rapsgelb-opak
B



36
gelb-fluoreszent
C



42
grün-opak
A B



45
lindgrün-opak
B



46
grün-fluoreszent
C



52
schwarz-opak
A B



62
weiß-opak
A B



72
beige-opak
A B



73
dunkelbraun-transparent
A



74
mittelbraun-transparent
A



101
türkis-metallic
D



102
mintgrün-metallic
D



103
lila-metallic
D



104
weinrot-metallic
D



105
silber-metallic
D



106
blau-metallic
D



107
limettengrün-metallic
D



108
rot-metallic
D



109
pink-metallic
D



110
gelb-metallic
D



111
orange-metallic
D



112
violett-metallic
D



113
schwarz-metallic
D